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# Instructions for the Use of Surgical Instruments

## ***Use of new Surgical Instruments***

Every instrument must be cleaned and sterilised before it is used.

## ***Caution***

The instruments may only be used for their intended purpose in the surgical specialities by educated and qualified personal. The surgeon shall be responsible for the proper selection of the instruments for each application, for obtaining the appropriate training, knowledge and experience and for their operative use.

Instruments4you as manufacturer and seller cannot accept any liability for immediate or consequential damages caused by inappropriate application and use or by inappropriate sterilisation and maintenance of the instruments.

If instruments are repaired by any companies or persons not authorized Instruments4you to do so, all warranties are becoming null and void.

Carefully examine each surgical instrument for breaks, cracks, deformations and malfunctions before use. It is especially essential to check areas such as blades, points, ends, stops and snaps as well as movable parts. Instruments that are worn out, corroded, deformed, porous or damaged in any other way, must be sorted out.

## ***Storage***

Instruments should be stored in a clean, dry, moisture free area. Instruments should be stored individually in their shipping carton or in a protective tray with partitions. Protect tips, edges etc. with tubing, protecting caps, gauze or fabric. Make sure that no chemicals are close to or in the storage area.

## ***Used materials***

Stainless steel as per DIN EN ISO 7153-1; Titanium alloys as per DIN ISO 5832-3; Pure titanium as per DIN ISO 5832-2; Light metal (Aluminium)

## ***Steel instruments***

The high-grade steels (rustproof, stainless) that are used for manufacturing surgical instruments create due to the chemical composition specific passive layers as protective surfaces. Those steels however are only to a certain extent resistant against attacks of chloride ions and aggressive waters! Chloride ions mainly cause pitting, but can also cause stress corrosion cracking. The greatest danger is water in which considerable quantities of common salts (sodium chloride) are dissolved.

## ***Titanium instruments***

Instruments made from pure titanium or titanium alloy can be handled and treated like steel instruments and no special precautions gave to be taken. Some titanium instruments are completely or partially anodised in blue colour for identification purpose.

## ***Aluminium instruments***

Only non-alkaline, neutral cleaning agents in combination with fully demineralised water must be used. Otherwise damages to the anodized surface are possible. Alkaline cleaning causes marks and colour fading on the surface particularly of coloured instruments already after just a few cycles.

In addition to the endeavours undertaken by the manufacturer with regards to the selection of the proper materials and its careful processing, the user has to ensure continuous and proper care of the surgical instruments as well as proper preparation, cleaning and sterilisation.

## ***We recommend the following methods and procedures for the preparation of our reusable surgical instruments:***

### ***Manual Cleaning***

The instruments must be disinfected and cleaned immediately after use. Contaminations on the instrument must not get dry or encrusted, as this could cause difficulties in cleaning and disinfection.

The following points are to be observed:

- Solutions used for the mechanical cleaning must be prepared strictly following the instructions given by the manufacturer
- For the cleaning of cannulas and dead-end holes a suitable brush must be used so any area can be reached.
- Remove blood and all contaminations with a soft brush and a mild neutral or alkaline (except for Aluminium!) detergent. Never use metal brushes or metal sponges for cleaning.
- To ensure proper function of the instrument, make sure that all movable parts are thoroughly cleaned.
- Clean instruments with hinges and box-locks in open as well as in closed position.
- Detach instruments for cleaning as far as possible.
- Pay special attention to the cleaning of slots, gaps, ratchets, box-locks, cannulations and dead-end holes.
- Surgical instruments should be placed in proper carriers, such as perforated trays, wire baskets etc.

### ***Ultrasonic treatment***

For ultrasonic treatment instruments should be placed in open condition on proper perforated trays or in wire baskets. Please ensure to avoid any "wave shadows" or covering surfaces caused by wire baskets or perforated trays or by large or bulky instruments.

Warm water without any additives does not have a satisfactory cleaning result and therefore a suitable cleaning agent should be added. Follow strictly the instructions given by the manufacturers regarding concentration.

The temperature of the detergents in the ultrasonic basin should range between 104°F and 113°F (40°C and 45°C). If the temperature is below 104°F (40°C) no cleansing effect is guaranteed. If the temperature is above 113°F (45°C), the proteins may coagulate.

A too dirty solution in the ultrasonic basin decreases the cleaning effect. Therefore, the solution should be renewed at intervals according to the instructions given by the manufacturer.

Ultrasonic waves for 3 minutes at a frequency of at least 35 kHz have proved to be sufficient in everyday use (Follow manufacturers' instructions.)

After ultrasonic treatment all instruments must be rinsed and checked for loose parts (e.g. screws etc.). For rinsing fully demineralised or distilled water must be used to avoid water spots.

# Instructions for the Use of Surgical Instruments

## **Chemical disinfection**

- The temperature of soaking solutions used for chemical disinfection should not exceed 140°F (60°C).
- Thinnings have to be made using fully demineralised or distilled pure water only. Detergent or cleaning agent must not be added. Follow precisely the instructions given by the manufacturer of the solution regarding dosage and induction time.
- The disinfection solutions must be refreshed daily. Reusing them can cause an increase of the dosage through evaporation (->corrosion risk) or a too high contamination level (->corrosion risk and reduced efficiency).
- After chemical disinfection all instruments must be rinsed with pure flowing water. To avoid water spots only fully demineralised or distilled water must be used.
- Dry surgical instruments immediately after each cleaning, disinfection and rinsing cycle.

## **Machine Cleaning and Disinfection**

- Machine cleaning and disinfection is always a preferable method compared to manual cleaning since machine procedures can be standardised.
- Follow the operating and loading instructions provided by the manufacturer of the washing machine. Use only the detergents and cleaning agents recommended by the manufacturer for the specific purpose.
- Hinged and box-lock instruments must be loaded and cleaned in open condition. Place instruments into the machine in a way that allows the water to flow out of cannulations and dead-end holes.
- Take instruments apart as much as possible for cleaning.
- Machine cleaning and disinfection is only suitable for instruments with long or thin cannulations if the hot disinfection solution can actually flow through them.
- When removing instruments from the washing machine, pay special attention to the proper cleaning of slots, gaps, ratchets, box-locks, cannulations and dead-end holes. Check for any visible remaining contaminations. If necessary clean manually and/or repeat cycle.

## **Steam sterilization / Autoclaving**

Sterilize all instruments before use.

Recommended sterilization method → Steam sterilization with fractionated vacuum according to EN 554  
Recommended temperature → 273°F (134°C); Recommended pressure → 3 bar; Leave on time → ≥ 5 min.

When using autoclaves for sterilization of surgical instruments, it has to be strictly ensured that the steam used is absolutely free of foreign substances such as corrosive particles or dirt to avoid subsequent corrosion or dirt (scum) deposit. Please observe strictly the instructions for use given by the manufacturers of autoclaves. Do not use any damaged instruments. Equipment and methods used for cleaning, disinfection and sterilization have to be in accordance with the following standards and recommendations:

DIN EN ISO 17664 → Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable devices

EN 285 → Sterilization – Steam-sterilizers – Large Sterilizers

EN 554 → Sterilization of medical devices – Validation and routine control of sterilization by moist heat

DIN EN 556-1 → Sterilization of medical devices – Requirements for medical devices to be designated “STERILE” – Part 1: Requirements for terminally sterilized medical devices

DIN 58946-7 → Sterilization, Steam-sterilizers, structural requirements for Large Sterilizers

Proper Maintenance of Instruments, Working group instrument preparation <http://www.a-k-i.org/>

Repeated instrument preparation has no significant influence on the lifetime of the surgical instrument. The lifetime of the instrument is usually determined by wear and tear or mutual damage during use.

## **Maintenance of instruments**

Maintenance of surgical instruments means lubrication with physiologically inert instrument oil (acc. to DAB 8 or Ph.Eur. or Usp) particularly of the joints. Make it a basic rule to thoroughly lubricate surgical instruments prior to checking for function. All movable parts (joints) and cutting blades of scissors have to be lubricated. This avoids metal abrasion when checking for function. Lubricants used must guarantee that even after frequently repeated use a “sticking” of joints through a multiplying effect is avoided.